



DispozABLE Spacer

Technical specifications

Tradition and Innovation - Visionary thinking and a fascination with technology have guided us throughout our history to develop innovative products of outstanding reliability. Anticipating trends to improve the quality of life.

DispozABLE Spacer

Technical specifications

DispozABLE™ Spacer

DispozABLE Spacer

| | |
|--------------------|---------------|
| Dimensions | ø84mm x 107mm |
| Weight | 15g |
| Volume | 227ml |
| Resistance to flow | Free flow |

Performance Data

Fine particle dose (salbutamol):

Delay 0 secs = 54.5mcg (SD 4.8)

Delay 1 secs = 37.7mcg (SD 7.6)

Product

Part Code

| | |
|-------------------------------|---------|
| DispozABLE Spacer (box of 10) | 3609550 |
|-------------------------------|---------|

Product information

DispozABLE Spacer has been designed to meet a specific market need. It is a low cost, fast response spacer for delivering medication promptly. It is designed to be used straight from the box without pre-washing in an emergency situation. For example in an ambulance or school. DispozABLE Spacer is to be used as a single use spacer. This simple intuitive design allows non-professionals to assemble it quickly. Clear instructions are included on the spacer, adding to the element of ease.

- Simple to use
- Effective
- For use in an emergency situation: ambulance or school
- Easy to assemble
- No mask required
- Can be used with tidal breathing technique, assisted by Flo-Tone.

DispozABLE Spacer is manufactured in accordance with all appropriate EU regulations.

© CLEMENT CLARKE INTERNATIONAL Ltd.

CLEMENT CLARKE INTERNATIONAL Ltd | Edinburgh Way | Harlow | Essex | CM20 2TT, GB
Phone +44 (0)1279 414969 | Fax +44 (0)1279 456300 | resp@clement-clarke.com | www.clement-clarke.com